


Novel Biomarkers of Cisplatin-induced Kidney Damage

D. Eisinger², J. Mapes², Salah-Dine Chibot¹, Gerard Maurer¹, F. Dieterle¹

¹Translational Sciences, Novartis, CH-4002 Basel, Switzerland; ²Rules-Based Medicine Inc. (RBM), Austin TX, USA



Methods

Control and Cisplatin study samples with previously documented proximal tubular injury were utilized to investigate novel biomarker candidates. Urine was collected at termination days 3, 7, 14 and 21, plasma at days 1, 3, 7, 21. Histopathology was performed with a standardized hierarchical lexicon. A panel of 58 protein biomarkers comprising cytokines, inflammatory markers, growth factors, and tissue factors were quantified in urine and plasma samples with the RBM RodentMAP[®] v2.0 multiplex immunoassay service. Performance of biomarker candidates from the panel was assessed with a Receiver Operating Characteristic (ROC) analysis.

RodentMAP[®] v2.0 | Antigens

1. Apolipoprotein A1	20. Interleukin-1 beta	39. MIP-2
2. C-Reactive Protein	21. Interleukin-2	40. MIP-3 beta
3. CD40	22. Interleukin-3	41. MMP-9
4. CD40 Ligand	23. Interleukin-4	42. MCP-1
5. Endothelin-1	24. Interleukin-5	43. MCP-3
6. Eotaxin	25. Interleukin-6	44. MCP-5
7. Epidermal Growth Factor	26. Interleukin-7	45. Myeloperoxidase
8. Factor VII	27. Interleukin-10	46. Myoglobin
9. Fibrinogen	28. Interleukin-11	47. Oncostatin M
10. FGF-basic	29. Interleukin-12p70	48. RANTES
11. FGF-9	30. Interleukin-17	49. Serum Amyloid P
12. GCP-2	31. KC/GRO alpha	50. SGOT
13. GM-CSF	32. Leukemia Inhibitory Factor	51. Stem Cell Factor
14. GST-alpha	33. Lymphotactin	52. Thrombopoietin
15. Haptoglobin	34. M-CSF	53. TIMP 1
16. Immunoglobulin A	35. MDC	54. Tissue Factor
17. Inducible Protein-10	36. MIP-1 alpha	55. Tumor Necrosis Factor-alpha
18. Interferon-gamma	37. MIP-1 beta	56. VCAM-1
19. Interleukin-1 alpha	38. MIP-1 gamma	57. VEGF
		58. von Willebrand Factor

Introduction

In 2008, the FDA and EMEA issued a joint framework allowing for the biomarker qualification and regulatory endorsement of 7 new urinary biomarkers of drug-induced kidney injury in rats. We sought to screen for other novel biomarkers of drug induced kidney injury by utilizing the original Novartis study samples that were the source for the Predictive Safety Testing Consortium's (PSTC) submission to the FDA and EMEA.

Results


Several new biomarkers not linked to previously described kidney biomarkers were significantly increased in urine and plasma and correlated to histopathological kidney damage severity. In addition, ROC analysis indicated that several of the novel analytes are sensitive and specific markers of drug-induced kidney damage with a significantly higher AUC of ROC than BUN and serum Creatinine.

ROC Analyses

Biomarker (Blinded)	Medium	Norm	Analysis Type	AUC	AUC Stdev	Dir	Thrd	Spec	Sens	N Ctrls	N Injur
Serum Creatinine	Serum		Exclusion	0.527	0.082	1	1.075	1.000	0.286	23	28
			Inclusion	0.586	0.065	1	1.065	0.955	0.276	67	29
BUN	Serum		Exclusion	0.615	0.079	-1	0.850	0.957	0.357	23	28
			Inclusion	0.520	0.065	-1	0.751	0.955	0.103	67	29
Protein J	Serum		Exclusion	0.845	0.055	1	1.764	0.957	0.643	23	28
			Inclusion	0.683	0.062	1	3.364	0.955	0.103	67	29
Protein K	Serum		Exclusion	0.828	0.058	1	1.149	0.957	0.536	23	28
			Inclusion	0.712	0.061	1	1.276	0.955	0.310	67	29
Protein L	Serum		Exclusion	0.781	0.064	1	1.592	0.957	0.571	23	28
			Inclusion	0.695	0.061	1	2.773	0.955	0.207	67	29
Protein M	Serum		Exclusion	0.762	0.066	1	1.417	0.957	0.571	23	28
			Inclusion	0.689	0.062	1	2.088	0.955	0.207	67	29
Protein N	Serum		Exclusion	0.792	0.063	1	1.397	0.957	0.5	23	28
			Inclusion	0.682	0.062	1	2.096	0.955	0.207	67	29
Protein O	Urine		Exclusion	0.846	0.059	1	1.413	1.000	0.538	17	26
			Inclusion	0.827	0.055	1	1.413	0.951	0.538	41	26
Protein P	Urine	UCrea	Exclusion	0.826	0.062	1	2.229	1.000	0.385	17	26
			Inclusion	0.845	0.053	1	1.524	0.952	0.654	42	26
Protein Q	Urine	UCrea	Exclusion	0.910	0.046	1	1.846	1.000	0.680	16	25
			Inclusion	0.922	0.039	1	1.656	0.974	0.720	39	25
Protein L	Urine	UCrea	Exclusion	0.889	0.050	1	1.574	1.000	0.808	17	26
			Inclusion	0.861	0.050	1	2.417	0.952	0.654	42	26
Protein M	Urine	UCrea	Exclusion	0.871	0.054	1	1.470	1.000	0.769	17	26
			Inclusion	0.845	0.053	1	2.281	0.952	0.538	42	26
Protein R	Urine	UCrea	Exclusion	0.844	0.059	1	1.755	1.000	0.538	17	26
			Inclusion	0.831	0.055	1	2.117	0.952	0.500	42	26
Protein S	Urine	UCrea	Exclusion	0.837	0.060	1	1.930	1.000	0.500	17	26
			Inclusion	0.855	0.051	1	1.738	0.952	0.577	42	26
Protein T	Urine	UCrea	Exclusion	0.819	0.064	-1	0.040	1.000	0.038	17	26
			Inclusion	0.796	0.059	-1	0.138	0.952	0.192	42	26
Protein U	Urine	UCrea	Exclusion	0.894	0.049	1	1.567	1.000	0.538	17	26
			Inclusion	0.877	0.048	1	1.588	0.952	0.500	42	26
Protein V	Urine		Exclusion	0.833	0.061	-1	0.144	1.000	0.154	17	26
			Inclusion	0.715	0.066	-1	0.286	0.952	0.154	42	26
Protein W	Urine		Exclusion	0.934	0.038	1	1.280	1.000	0.808	17	26
			Inclusion	0.934	0.036	1	1.280	0.952	0.808	42	26
Protein X	Urine		Exclusion	0.846	0.059	1	1.356	1.000	0.462	17	26
			Inclusion	0.810	0.057	1	1.423	0.951	0.346	41	26
Protein Y	Urine		Exclusion	0.887	0.050	1	1.183	1.000	0.577	17	26
			Inclusion	0.870	0.049	1	1.144	0.952	0.577	42	26

Plots

The plots show the biomarker levels for all samples investigated. All values are represented as fold-changes. The animals are ordered by dose-group (with increasing doses) and within each dose-group by termination time point (with increasing time). The symbols and the colors represent the histopathology readout for proximal tubular damage for the corresponding animal (red = no histopathology finding observed, green = grade 1, black = grade 2 on a 5 grade scale). The data point labels represent the termination time point and the x-axis labels represent the dose level.



Conclusions

The newly identified biomarkers have potential to complement the already known and partially qualified biomarkers for drug induced kidney injury in particular for monitoring inflammatory processes, whereas the 7 endorsed biomarkers are mechanistically linked to functional changes, cell injury and regeneration processes.